



Measurement of Volatile Organic Compounds in Indoor Air of Museums by Serially Connected Passive Samplers

M. Butsugan¹ and Y. Sekine²

¹Hitachi Chemical Co. Ltd

²School of Science, Tokai University




Passive samplers

- Simple device suitable for assessing indoor air quality and personal exposure
 - Previous passive samplers;
 - DSD-DNPH
for aldehydes and ketones
 - VOC-SD
for VOCs
- (available from SIGMA-ALDRICH, Japan)






DSD-DNPH



Schematic view of the sampler

DNPH: 2,4-dinitrophenylhydrazine




SEM image of the diffusion filter

Porosity: 34.5%

Air cavities: ~20 μ m

Determination by HPLC




Reaction of carbonyl compounds with DNPH

DNPH derivatives are extracted by 10ml of acetonitrile and subsequently determined by reverse-phase HPLC with UV-Vis detector at 315nm.





VOC-SD



Schematic view of the sampler

Trapped VOCs are extracted by 2ml of carbon disulfide and subsequently determined by GC/FID or GC/MS.





Sampling rates

- Sampling rates of both passive samplers were practically determined by small chamber experiments comparing with active samplings

Sampler	Analyte	Sampling rate, α (ml/min)	ref.
DSD-DNPH	formaldehyde	71.9	(1)
	acetaldehyde	59.4	(1)
	acetone	51.7	(1)
	glutaraldehyde	40	(2)
VOC-SD	benzene	47	(3)
	toluene	42	(3)
	ethylebenzene	32	(3)
	<i>m-,p</i> -xylene	35	(3)
	styrene	37	(3)
	<i>o</i> -xylene	32	(3)

(1) Uchiyama et al., *Atmos. Environ.*,38,6319(2004) (2) Sekine et al.,*J. Health Sci.*,51,629(2005)

(3)Butsugan et al., *Proc.of Indoor Air 2005*,2289(2005)

Serially connected sampler



- The sampling portions of DSD-DNPH and VOC-SD are connected in series to give a single tube (DSD-Serial) for samplings of multiple indoor air pollutants.
- The new sampler was applied to field measurements of indoor air pollutants in a museum and validated the sampling performance.


● ● ● | Field tests in museum

- Takehiko Miyanga Memorial Museum, Hadano, Kanagawa, Japan
- Major displays: Oil paintings
- Floor area : 192m²
(exhibition room)
- 24hrs-air ventilation system
Air exchange rate:>0.7 ACH
- Sampling period: July 9 –11, 2006 (Museum was closed.)
- Sampling duration: 38.5h
- Temp: 24.5°C, R.H.70%





Sampling sites




Number of samplers deployed at each sampling site

Sampler	A	B	C	D	travel blank
DSD-Serial	6	2	2	2	3
DSD-DNPH	6	-	-	-	3



Deployment of the DSD-Serial



Deployed at the height of 2 m from the floor



Analytical procedure

<Aldehydes & ketones>

DNPH coated silica gel



Transfer to vial



← 3ml of acetonitrile

Extraction with mild shaking



Centrifugation (3000rpm,5min)



HPLC analysis
(Hitachi Elite System)

<VOCs>

Carbon molecular sieves



Transfer to vial



← 2ml of carbon disulfide

Extraction at 5°C for 2hrs



Centrifugation (3000rpm,5min)



GC/MS analysis
(Shimadzu QP-5000 System)

Comparison of collection amount of aldehydes and ketones




Figure. Comparison of collection amount of formaldehyde and acetone by co-located two kinds of samplers at site A (n=6).

No significant difference was found (t -test, $p < 0.05$)




Travel blanks and LOQ

Analyte	Travel blank (μg)	LOQ(ppb)
formaldehyde	0.072 \pm 0.004	0.14 *
acetaldehyde	0.063 \pm 0.006	0.20 *
acetone	0.17 \pm 0.006	0.17 *
benzene	not detected	0.69 **
toluene	0.13 \pm 0.05831	1.6 *
<i>m</i> -, <i>p</i> -xylene	0.062 \pm 0.01301	0.37 *
ethylebenzene	not detected	0.71 **
styrene	not detected	0.66 **
<i>o</i> -xylene	not detected	0.57 **
<i>p</i> -dichlorobenzene	not detected	0.42 **


LOQ was defined as, * : 10σ , ** : S/N=10

Sampling duration: 38.5h



Formaldehyde in the museum


formaldehyde (IAQ guideline=80ppb)



Possible emission source: wooden materials for flooring etc.



Acetaldehyde in the museum




Possible emission source: wooden materials for flooring etc.



Toluene in the museum

toluene(IAQ guideline=70ppb)




Possible emission source: background outdoor air



m,p-xylene in the museum

m,p-xylene (IAQ guideline=200ppb)




Possible emission source: oil paintings



p-dichlorobenzene in the museum

p-dichlorobenzene (IAQ guideline=40ppb)



Possible emission source: background outdoor air



Conclusion

We have developed a serially connected passive sampler (DSD-Serial) and applied to field measurements of indoor air pollutants in the museum and lead the conclusions as follows.

- Sampling rates of previous passive samplers (DSD-DNPH & VOC-SD) can be used for the DSD-Serial, because no significant difference was found in the collection amounts of analytes.
- DSD-Serial showed a high sensitivity enough to assess indoor air quality of a museum.
- Indoor air quality of the Takehiko Miyanaga Memorial Museum was well controlled by mechanical ventilations.



Acknowledgement

Mr Hitoshi Kimura

Hadano City Government Board of Education
Miyanaga Takehiko Memorial Museum

Mr Daisuke Oikawa

Preventive Clinical Laboratories, Ltd.

Ms Yukiko Kado

School of Science, Tokai University

Mr Toshimitsu Kuriwada

Project manager, Tokai University